

PRICE SCHEDULE /QUOTATION

DNA Testing Option	DNA test for personal Knowledge	DNA test for Legal purpose	Lab process time
Paternity Test. Mother, Child and Alleged Father OR Child and Alleged Father without mother	20,000 Ksh	20,000 Ksh	5-10 Working days
Maternity Test. Mother and child	20,000 Ksh	20,000 Ksh	5-10 Working days
Relationship Testing			
Full Siblingship. Two alleged full siblings with or without common parent	20,000 Ksh	20,000 Ksh	5-10 working days
Half sibling ship. Two alleged half siblings with or without common parents	20,000 Ksh	20,000 Ksh	5-10 Working days
DNA Profile Individual	10,000	10,000	5-10 working days
DNA Forensic Viability Test			
A test on unusual samples is used to detect presence of usable DNA. If the viability test reveals the presence of usable DNA, A standard DNA test will then proceed.	Fee is variable and depends on sample type and condition.		Duration depends on type of sample
Child ID Services Nursery Confirmation. 1 child, legally binding test	20,000	20,000	5-10 working days

*Additional tested party Ksh 10,000

VIABILITY TESTING

Unusual samples include

Percentage (%) means rate of success

Category A \geq 90%	Category B 60–90%	Category C 20–60%	Category D \leq 20%
10,000 Ksh	10, 000 Ksh	12,500 Ksh	12,500 Ksh
Blood fresh or in cold storage less than 30 days old	Blood / Semen stain on clothing includes bandages, gauze, and facial tissue paper or napkins	Drinking Straw	Hair with roots ensure roots are present)7-10 pieces) shed hair will not work
Blood on filter paper as called whatman, Matrix or parchment less than 30 days old. Blood on filter paper as called whatman, Matrix or parchment less than 30 days old.	Tissue (muscles, organs) depends on age & storage conditions, non-embalmed	Feminine pad	Tooth Brush
Blood on filter paper as called whatman, Matrix or parchment less than 30 days old.	Dried umbilical cords must have a reference sample from mother.	Soda can / drinking glass swab well with cotton swab	Cigarette butts
Semen fresh or frozen	Cotton Swab (ear swab)	Plastic fork / spoon Swab well	Condom depends on storage condition / swab inside of condom